

Cascades of security

BROEN-LAB Emergency Shower Systems

General instructions for emergency showers
To be filled out by the installer and retained by the customer

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BROEN-LAB

Emergency Shower

BROEN-LAB emergency showers are used as first aid in accidents causing burns, scalds, chemical contamination etc. on the human body or eye.

These emergency showers comply with EN 15154 and ANSI Z 358.1 norms.

The wording "emergency shower" refers in the following to the complete range of body-, eye- and hand showers!

Placement of Emergency Shower

The BROEN-LAB emergency shower must be placed as close to the person/place of work as possible. Preferably it must be placed so that an injured person is able to reach the shower within 10 seconds corresponding to a max. distance of 10 to 15 metres (33 to 50 ft.). In case of strong acids or caustics the shower must be immediately adjacent to the hazard.

The shower must be easy to find, and right below the shower head there must be a free space of 0.8 m (32") in diameter. Marking on the floor of this zone is advisable in order to prevent the placing of goods, items or similar below the shower. Install the emergency shower close to a floor outlet.

Water Supply

Mains water only! Drinking water or water of equivalent quality.

Min. dynamic pressure: as indicated on the sticker on the product (must be obtained during operation)

Max. stat. and dyn. pressure: 10 bar (145 psi).

Min. flow rate at min. dyn. pressure:

Body showers:

EN 15154-1: 60 l/min

ANSI Z 358.1: 76 l/min (20 US gallons/min).

Eye and eye/face wash: 6 l/min per outlet (1.6 US gallons/min.)
i.e. min. 12 l/min (3.2 US gallons/min.) if 2 heads.

Ideal temperature:

Between 15°C / 59°F and 25°C / 77°F.

Please note i.e. ANSI Z 358.1: Tepid water!

The above indications are to be tested by measuring the pressure at the shower inlet and by measuring the flow rate. Water temperature is measured with a temperature gauge.

BROEN-LAB recommend having sufficient water for 30 minutes flushing.

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
BROEN-LAB anbefaler, at 30 minutters kontinuerlig skylning kan gennemføres.

Mounting

- **IMPORTANT - Flush the pipeline system to remove the impurities**
- Mount the BROEN-LAB emergency shower. See next paragraph for rec. dimensions etc.
- Please refer to the below illustration when mounting the BALLOFIX ball valve. It is important that unauthorized shut off is prevented.
- Connect the shower to the water supply

- Test the BROEN-LAB emergency shower.

Mounting recommendation

ANSI Z 358.1

- The top of the water jet of the activated eye shower must be 83,8-134,6 cm (33-53 in) above ground and minimum 15,3 cm (6 in) from the wall or other obstructions
- The pull rod must be max. 175,3 cm (69 in) above ground
- Bottom tip of shower head must be 208,3-243,8 cm (82-96 in) above ground

EN 15154 part 1

- Bottom tip of the showerhead on a body shower must be 220 +/- 10 cm above floor


Maintenance

At least every month the operation of the shower should be checked. Beyond this test the BROEN-LAB emergency showers are maintenance-free.

Acc. to ANSI self-contained emergency showers shall be visually checked weekly to determine if the water that flushes the shower needs to be changed/supplemented

Maintenance of Filter BALLOFIX[®] if mounted

- The ball is turned clockwise to its closed position.
- The side plug is removed with a coin or a screwdriver.
- The filter is removed by a pair of tongs.


- The filter can be cleaned with compressed air.

If you have any problems or questions regarding BROEN-LAB emergency showers, please get in touch with your local BROEN-LAB distributor.

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
GB

Signs

Included in the packaging of all BROEN-LAB emergency showers the customer receives an international emergency symbol.

This symbol must be placed where it is visible and where it clearly indicates the installation of an emergency shower.

Instructions for Use

- Start flushing immediately – the seconds counts
- Remove contaminated clothing which is not stuck to the skin
- Continue flushing
- In case of an eye injury keep the eyelid open when continuing flushing
- Call for an ambulance. Ask the doctor for advise
- Continue flushing in the ambulance
- Keep the flushing in the emergency room and during transport to the special ward

Do remember to inform all employees about:

- How to prevent accidents/injuries
- How to react should an emergency arise
- The location of the emergency shower
- How to operate the emergency shower
- How to help injured persons in the best way.

Warning!

Before the installation of your BROEN-LAB emergency shower, all debris which may cause leaks must be thoroughly flushed from the system. Should a fault appear in this product, and the system was found not to have been flushed out, the warranty will be void.

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