


SUSPENDED FITTING



SPECIFICATIONS




Pressure regulator with integrated shut-off valve for suspended mounting.

Inlet connection: female 1/4" NPT
(see ordering information for other options)


Outlet connection: female 1/4" NPT
(see ordering information for other options)

MOUNTING INFORMATION

Installation example



Drilling template



ORDERING INFORMATION

Brass: 2503... - models
Stainless Steel: 2504... -models

Fittings

Outlet pressure:

	without flow valve	with flow valve
0-1 bar	2503 710 10xx-yy	2503 711 10xx-yy
0-1.5 bar	2503 720 10xx-yy	2503 721 10xx-yy
0-2.5 bar	2503 730 10xx-yy	2503 731 10xx-yy
0-5 bar	2503 740 10xx-yy	2503 741 10xx-yy
0-10 bar	2503 750 10xx-yy	2503 751 10xx-yy
0-16 bar	2503 760 10xx-yy	2503 761 10xx-yy
0-25 bar	2503 770 10xx-yy	2503 771 10xx-yy
0-65 bar	2503 780 10xx-yy	2503 781 10xx-yy

Delivered with mounting accessories

Inlet/outlet options, tube fittings (to be ordered separately)

	Straight Brass	Straight Stainless Steel	Elbow Brass	Elbow Stainless Steel
OD 1/8" x male 1/4" NPT	96i01400	96i01040	96i01560	96i01100
OD 1/4" x male 1/4" NPT	96i00390	96i01030	96i01550	96i01090
OD 6mm x male 1/4" NPT	96i01120	96i01020	96i01320	96i01080
OD 8mm x male 1/4" NPT	96i00380	96i01110	96i01540	96i00430
OD 10mm x male 1/4" NPT	96i01240	96i01430	96i01390	96i01680

For further information please contact your local BROEN-LAB representative.